







DiaDent Group International www.diadent.com | 1.877.342.3368 Follow us on FACEBOOK





# THE X FACTORS:

#### **Enhanced Safety**

- Convex triangular x-section reduces rotational fiction and contact with canal wall
- Progressive decreasing % taper design enhances flexibility and safety
- Navigate calcified or severely curved canals with confidence and control
- Fracture resistance is 143% greater than Dia-PT™ Files
- High corrosion resistance with coated surface

### Greater accessibility


- Heat treated premium nickel titanium construction offers extreme flexibility
- 11 mm length handle offers improved accessibility to teeth
- Flat non-cutting file tip efficiently removes debris and soft tissue

#### **Highly Integrable**

- Designed to work interchangeably with your current technique and system
- Predictable results every time
- Perfect shape-matching gutta percha and paper points: Dia-ProT™

#### ORDER INFORMATION

Size	19MM	21MM	25MM
#DX	571-100		
#D1		571-101	571-201
#D2		571-102	571-202
#D3		571-103	571-203
#D4		571-104	571-204
#D5		571-105	571-205
#Asst D1/D4		571-191	571-291
			(Pkg. of 4)






# Dia-X<sup>TM</sup> File

•RECOMMENDED SPEED: 250~350 RPM

•RECOMMENDED TORQUE: 2~3 NCM





2



Before the endo procedure, disinfect inside of the canals using EDTA. Use DX File to shape the coronal section and enlarge the opening of the canal.

## K-File




Determine the working length with an X-ray and scout the canal space using K-Files #8 and #10.

### K-File



Follow the glide path using K-File #15 and enlarge the opening of the canal.


4



D1 expands the glide path and enlarges to widen the coronal section. Irrigate with NaOCI before engaging the file. Shape the middle part of the coronal section using D1 File.



Continue shaping the middle and apical portions with D1 file using gentle force to brush laterally. Cut the dentin on the outstroke to improve straight-line access and apical progression.




After providing adequate irrigation, continue shaping the rest of the section (up to the apex) using D1 File.


7




D2 enlarges to the apical 2/3. Shape the canal section with a gentle brushing action using D2 File. Confirm working length frequently, especially when shaping severely curved canals.



Shape down to the apex using D<sub>3</sub> File with an in-and-out pecking motion. Finish shaping when it matches the size of K-File #20.



If the apex is larger than hand file #20, continue shaping with D4 (if #20 hand file is loose at length) and/or D5 (if #25 hand file is loose at length).





Dia-X DX: Access File (19mm) Dia-X D1: ProTaper® S1

> Dia-X D2: ProTaper® S2 Dia-X D3: ProTaper® F1

> Dia-X D4: ProTaper® F2

Dia-X D5: ProTaper® F3

Use exactly like ProTaper® Gold

### Dia-ProT™ Gutta Percha & Paper Points (Millimeter-Marked)





### **TIPS**

- Inspect files for distortion or wear, such as uneven flutes, dull spots or bending
- Lubricate files generously with Dia-Prep  $^{\text{TM}}$  Pro
- Remove the debris from the flutes before continuing
- Irrigate and re-irrigate after using each rotary file
- Use very light pressure. Never force instruments to working length
- Use the files for no more than one or two seconds in the canal
- Autoclave files after each use

### Dia-ProT Gutta Percha (60/box)

ML150-S601 F1 (matches D3) ML<sub>150</sub>-S602 F<sub>2</sub> (matches D<sub>4</sub>)

ML150-S603 F3 (matches D5)

ML150-S691 Asst (F1/F2/F3)

### **Dia-ProT Paper Points**

MP250-S601 F1

MP250-S602 F2

(100/box)

MP250-S603 F3

MP250-S691 Asst(F1/F2/F3)