

# DIA-PROISO .04/.06 PLUS

## MILLIMETER MARKED SPECIAL TAPER PAPER POINTS

- Millimeter marked Special Taper Absorbency Paper Points.
- Designed to be used for drying and cleaning the canals prepared with **Profile, Hero & K3 files.**
- Available in .04 & .06 taper.




## Available Sizes

Size	Dia-ProISO .04 Plus	Dia-ProISO .06 Plus
	P-63	P-63
#15	MP 258-603	MP 259-603
#20	MP 258-604	MP 259-604
#25	MP 258-605	MP 259-605
#30	MP 258-606	MP 259-606
#35	MP 258-607	MP 259-607
#40	MP 258-608	MP 259-608
#45	MP 258-609	MP 259-609
#50	MP 258-610	MP 259-610
#55	MP 258-611	MP 259-611
#60	MP 258-612	MP 259-612
#70	MP 258-613	MP 259-613
#80	MP 258-614	MP 259-614
#15/40	MP 258-691	MP 259-691
#45/80	MP 258-692	MP 259-692

(Pkg. of 100)

## DIA-PROISO .04/.06 PLUS IN THE CELL PACKAGE




- Points remain sterile until opened

Cell pkg	Dia-ProISO .04 Plus	Dia-ProISO .06 Plus
#15	MP258-703	MP259-703
#20	MP258-704	MP259-704
#25	MP258-705	MP259-705
#30	MP258-706	MP259-706
#35	MP258-707	MP259-707
#40	MP258-708	MP259-708
#15/40	MP258-791	MP259-791

(Pkg. of 120)

## DIA-PROT



- Designed to be used to complement **Protaper Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#F1	MP250-S601
#F2	MP250-S602
#F3	MP250-S603
#F4	MP250-S604
#F5	MP250-S605
#F1/F3	MP250-S691
#F4/F5	MP250-S692
#F1/F5	MP250-S693

(Pkg. of 100)

## DIA-PROT NEXT\* NEW




- Designed to be used to complement **Protaper Next Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#X2	MP267-S602
#X3	MP267-S603
#X4	MP267-S604
#X5	MP267-S605
#X2/X5	MP267-S691
#X2/X3	MP267-S692

(Pkg. of 100)

## DIA-PRO R




- Designed to be used to complement **RECIPROC® Files**.
- Millimeter Marked Points

Size	ITEM No.
#25	MP264-605
#40	MP264-608
#50	MP264-610
#Asst(25/40/50)	MP264-691

(Pkg. of 100)

## DIA-PRO W



- Designed to be used to complement **Wave One Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	MP265-S601
#Primary	MP265-S602
#Large	MP265-S603
#Asst(Small/Primary/Large)	MP265-S691

(Pkg. of 100)

## DIA-PRO W GOLD\* NEW



- Designed to be used to complement **Wave One Gold Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	MP268-S601
#Primary	MP268-S602
#Medium	MP268-S603
#Large	MP268-S604
#Asst(Small/Primary/Medium/Large)	MP268-S691

(Pkg. of 100)

## PACKAGE INFORMATION

### • P-63 Spill Proof Package

6 compartments with outer paper box.  
Box dimension: 73 x 17 x 49mm  
(W x H x D)



### • P-43 Vial Package

Contains 6 light weight plastic square vials with color coded caps  
Box dimension: 77 x 17 x 50mm  
(W x H x D)



### • SP-60 Slide Package

One hand slide open lid  
Cellophane shrink-wrapped package  
Box dimension: 70 x 13 x 45 mm  
(W x H x D)



### • P-33 Cell Package

Points remain sterile until opened  
Box dimension: 78 x 19 x 50mm  
(W x H x D)

